



PASS BOX BSPB-101

PASS BOX BSPB-101

It is used to pass goods between room of different level of cleanliness , so as to avoid air convention between them and minimize pollution extent

Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.

BSPB-101 PASS BOX



Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness

SPECIFICATIONS

Model	BSPB-101
UV Lamp	8Wx1
External Size (WxDxH)	660x570x580 mm
Internal Size (WxDxH)	500x500x500 mm
Material	304 stainless steel
Interlock	Electronical
Gross Weight	55 kg
Packing Size(WxDxH)	780x700x770 mm
Power Supply	220V, 50Hz; 110V, 50/60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com