



RECIPROCAL SHAKER BSRP-102

RECIPROCAL SHAKER BSRP-102

Light weight and portable reciprocal shaker with mild to vigorous shaking suitable for a variety of applications. Integral electronic timer handles both timed and continuous operations. Low noise and stable operation with over speed protection make it an ideal choice in different labs. Interchangeable platforms enable the use of different vessels. Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory.

Also known as Laboratory Shaker.

BSRP-102 RECIPROCAL SHAKER

Compact space-saving design fits easily in small spaces

Brushless DC motor ensures maintenance-free, reliable and quiet operation

Electronic time switching clock controls time or continuous operation mode

Simultaneous display of speed, timer and operating mode

RS232 interface enables data logging & computer control

Over speed detection and protection

Accommodate a variety of platforms, non-slip rubber mat, conical flasks and other sample containers

SPECIFICATIONS

Model	BSRP-102
Speed Range	0-200 rpm
Temperature	5-40°C
Shaking Orbit	20 mm
Motion	Reciprocal
Maximum Load Capacity	3 kg
Operation	Continuous
Overall Dimension	330x270x110 mm
Display	Analog
Weight	3.1 kg
Power	30 W
Power Supply	240 V, 50/60 Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com