



## UVC STERILIZER BSTE-301

## UVC STERILIZER BSTE-301

UVC Sterilizer is cost effective, Chemical free, efficient and environment friendly used to kill bacteria, viruses, mold, fungi and micro-organisms.

Used in Clinics, Disinfection of Microorganisms, Research, Pharmaceuticals.

Also known as Germicidal Lamp.

## BSTE-301 UVC STERILIZER



Simple, easy and convenient-Connected to a wireless network for operation via a mobile phone or tablet devices

Compact and lightweight design-Easy to setup and flexible for any space

Effective disinfection-By using UVC germicidal lamps but no chemicals, it can effectively kill 99.9% germs, also no need to clean up after use

Schedule and timer-Easy to use a mobile phone or tablet to set up timer or schedule for disinfection

## SPECIFICATIONS

|   |                                  |
|---|----------------------------------|
| Model   | BSTE-301                         |
| Wavelength  | 253.7 nm                         |
| Lamp Type   | UVC Lamp                         |
| Intensity   | 1360 uW/cm <sup>2</sup>          |
| Lamp Holder   | 2G11 Base                        |
| Wi-Fi   | IEEE 802.11 b/g/n 2.4GHz         |
| Android 4.1, iOS 9.0 or higher version  | OS Language                      |
| English, French, Italian, German, Bulgarian, Czech, Danish, Greek, Polish, Portuguese, Russian, Slovak, Spanish, Traditional Chinese, Turkish | Dimension (WxHxD) and Weight(kg) |
| 260x170x475 mm 2.9 Kg   | Package Dimension and Weight     |
| 495x235x165 mm 3.3 Kg   | Power                            |
| 144 Watt  | Input                            |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)