



UVC STERILIZER BSTE-302

UVC STERILIZER BSTE-302

UVC Sterilizer is cost effective, Chemical free, efficient and environment friendly used to kill bacteria, viruses, mold, fungi and micro-organisms.

Used in Clinics, Disinfection of Microorganisms, Research, Pharmaceuticals.

Also known as Germicidal Lamp.

BSTE-302 UVC STERILIZER



Simple, easy and convenient-Connected to a wireless network for operation via a mobile phone or tablet devices

Compact and lightweight design-Easy to setup and flexible for any space

Effective disinfection-By using UVC germicidal lamps but no chemicals, it can effectively kill 99.9% germs, also no need to clean up after use

Schedule and timer-Easy to use a mobile phone or tablet to set up timer or schedule for disinfection

SPECIFICATIONS

| | |
|---|----------------------------------|
| Model | BSTE-302 |
| Wavelength | 253.7 nm |
| Lamp Type | UVC Lamp |
| Intensity | 620 uW/cm ² |
| Lamp Holder | 2G11 Base |
| Wi-Fi | IEEE 802.11 b/g/n 2.4GHz |
| Android 4.1, iOS 9.0 or higher version | OS Language |
| English, French, Italian, German, Bulgarian, Czech, Danish, Greek, Polish, Portuguese, Russian, Slovak, Spanish, Traditional Chinese, Turkish | Dimension (WxHxD) and Weight(kg) |
| 260x85x475 mm 1.8 Kg | Package Dimension and Weight |
| 470x209x138 mm 2.4 Kg | Power |
| 72 Watt | Input |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com