




## BIOSAFETY TRANSPORT BOX BTBT-202

## BIOSAFETY TRANSPORT BOX BTBT-202

The Biosafety Transport Box is a non-electric cooler used for the transportation of biological products requiring physical insulation, such as biological samples, pathogenic microorganism (virus) species, blood, etc.

Used in It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism (virus) species, blood, etc.

## BTBT-202 BIOSAFETY TRANSPORT BOX



## SPECIFICATIONS

|                    |   |
|--------------------|---|
| Model              | BTBT-202                                    |
| Capacity           | 80 L  |
| Temperature Range  | -25 to 30°C                                 |
| Alarm              | High&low temp, control error, power failure |
| External Dimension | 930x580x570 mm                              |
| Power              | DC12V; AC100~220V, 50/60Hz                  |
| Optional Parts     | Basket; Datalogger                          |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)