



## BIOSAFETY TRANSPORT BOX BTBT-202

## BIOSAFETY TRANSPORT BOX BTBT-202

The Biosafety Transport Box is a non-electric cooler used for the transportation of biological products requiring physical insulation, such as biological samples, pathogenic microorganism (virus) species, blood, etc.

Used in It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism (virus) species, blood, etc.

## BTBT-202 BIOSAFETY TRANSPORT BOX



## SPECIFICATIONS

Model	BTBT-202
Capacity	80 L
Temperature Range	-25 to 30°C
Alarm	High&low temp, control error, power failure
External Dimension	930x580x570 mm
Power	DC12V; AC100~220V, 50/60Hz
Optional Parts	Basket; Datalogger



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)