



DIGITAL DESKTOP ULTRASONIC CLEANER BULC-101

DIGITAL DESKTOP ULTRASONIC CLEANER BULC-101

Digital ultrasonic cleaner is effective way to clean lab tools. It is a high efficiency cleaner specifically developed to clean surgical instruments.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-101 DIGITAL DESKTOP ULTRASONIC CLEANER



High strength and environmental ABS injection molding
Top quality ultrasonic ceramics parts, quality higher than competitors
Industrial IC, high voltage isolation design, safe using and long working life
Damp proof and anti corrosion process of the PCB
Strict production and quality test control
Digital control, multi stage ultrasonic cleaning adjustment
The tank is polished by SUS304 stainless steel stamping point
2500V high voltage test to ensure the product safety, strict production process, aging process and quality control

SPECIFICATIONS

Model	BULC-101
Capacity	600 ml
Working frequency	42 Khz
Ultrasonic Power	35 W
Heating power	No
Timing time	1-5 mins
Temperature range	No
Cover	Yes
Drain valve	No
Washing basket	Standard
Inner tank dimension	155x95x48 mm
Packing size	247x190x145 mm
Gross weight	1.0 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com