



## DIGITAL DESKTOP ULTRASONIC CLEANER BULC-102

## DIGITAL DESKTOP ULTRASONIC CLEANER BULC-102

Digital ultrasonic cleaner is effective way to clean lab tools. It is a high efficiency cleaner specifically developed to clean surgical instruments.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

### BULC-102 DIGITAL DESKTOP ULTRASONIC CLEANER



- High strength and environmental ABS injection molding
- Top quality ultrasonic ceramics parts, quality higher than competitors
- Industrial IC, high voltage isolation design, safe using and long working life
- Damp proof and anti corrosion process of the PCB
- Strict production and quality test control
- Digital control, multi stage ultrasonic cleaning adjustment
- The tank is polished by SUS304 stainless steel stamping point
- 2500V high voltage test to ensure the product safety, strict production process, aging process and quality control

## SPECIFICATIONS

Model	BULC-102
Capacity	750 ml
Working frequency	42 Khz
Ultrasonic Power	35 W
Heating power	No
Timing time	1-5 mins
Temperature range	No
Cover	Yes
Drain valve	No
Washing basket	Standard
Inner tank dimension	150x126x44 mm
Packing size	253x224x150 mm
Gross weight	1.17 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)