



## DIGITAL DESKTOP ULTRASONIC CLEANER BULC-103

## DIGITAL DESKTOP ULTRASONIC CLEANER BULC-103

Digital ultrasonic cleaner is effective way to clean lab tools. It is a high efficiency cleaner specifically developed to clean surgical instruments.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

### BULC-103 DIGITAL DESKTOP ULTRASONIC CLEANER



High strength and environmental ABS injection molding  
Top quality ultrasonic ceramics parts, quality higher than competitors  
Industrial IC, high voltage isolation design, safe using and long working life  
Damp proof and anti corrosion process of the PCB  
Strict production and quality test control  
Digital control, multi stage ultrasonic cleaning adjustment  
The tank is polished by SUS304 stainless steel stamping point  
2500V high voltage test to ensure the product safety, strict production process, aging process and quality control

## SPECIFICATIONS

Model	BULC-103
Capacity	1400 ml
Working frequency	40 Khz
Ultrasonic Power	60 W
Heating power	No
Timing time	1-5 mins
Temperature range	No
Cover	Yes
Drain valve	No
Washing basket	Standard
Inner tank dimension	168x136x62 mm
Packing size	316x313x213 mm
Gross weight	2.25 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)