



# DIGITAL DESKTOP ULTRASONIC CLEANER BULC-104

## DIGITAL DESKTOP ULTRASONIC CLEANER BULC-104

Digital ultrasonic cleaner is effective way to clean lab tools. It is a high efficiency cleaner specifically developed to clean surgical instruments.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

### BULC-104 DIGITAL DESKTOP ULTRASONIC CLEANER

High strength and environmental ABS injection molding  
Top quality ultrasonic ceramics parts, quality higher than competitors  
Industrial IC, high voltage isolation design, safe using and long working life  
Damp proof and anti corrosion process of the PCB  
Strict production and quality test control  
Digital control, multi stage ultrasonic cleaning adjustment  
The tank is polished by SUS304 stainless steel stamping point  
2500V high voltage test to ensure the product safety, strict production process, aging process and quality control

## SPECIFICATIONS

Model	BULC-104
Capacity	2500 ml
Working frequency	40 Khz
Ultrasonic Power	70 W
Heating power	100 W
Timing time	1-8 mins
Temperature range	65°C
Cover	Yes
Drain valve	No
Washing basket	Standard
Inner tank dimension	247x150x78 mm
Packing size	325x245x210 mm
Gross weight	3.04 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)