



# DESKTOP MULTI FREQUENCY TYPE ULTRASONIC CLEANER BULC-201

# DESKTOP MULTI FREQUENCY TYPE ULTRASONIC CLEANER BULC-201

Desktop Multi frequency type Ultrasonic Cleaner uses multi frequency of 40/80/120 Khz for carrying out vigorous cleaning to remove dirt and heavy contaminants.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-201 DESKTOP MULTI FREQUENCY TYPE ULTRASONIC CLEANER



Digital display and control

Special multi frequency ultrasonic transducer, cleaning conversion efficiency is higher

No stud welding technology, cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

Whole machine noise reduction processing

All stainless steel material inner tank is made by punching SUS304 1.0mm

Multi frequency or three frequency cleaning modes can be automatically/ manually switched

2500V high voltage test to ensure the product safety, strict production, aging production quality control

## SPECIFICATIONS

Model	BULC-201
Capacity	11 L
Working frequency	28/40 Khz or 40/80 Khz or 40/80/120 Khz
Ultrasonic Power	240 W
Heating power	300 W
Timing time	1-30 mins
Temperature range	0- 80°C
Cover	Yes
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	300x240x150 mm
Packing size	412x357x380 mm
Gross weight	9.7 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)