



HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER BULC-503

HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER BULC-503

High frequency ultrasonic cleaner is mainly used in medical industry, computers, microelectronics computers. It uses high frequency to clean the sub micron dirt without damaging the equipment

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-503 HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER



Digital display and control ultrasonic power adjustable range from 0 to 100%

Special high frequency ultrasonic transducer (53/68/80/100/132 KHz)

No stud- welding technology, Cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

All stainless steel material inner tank is made by punching SUS304 1.0mm

2500V high voltage test to ensure the product safety, strict production, aging production quality control

SPECIFICATIONS

Model	BULC-503
Capacity	4.5 L
Working frequency	53/68/80 KHz
Ultrasonic Power	180 W
Heating power	200 W
Timing time	1- 99 mins
Temperature range	0- 80°C
Power regulation	No
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	300x150x100 mm
Packing size	300x150x100 mm
Gross weight	7.43 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com