



## HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER BULC-506

# HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER BULC-506

High frequency ultrasonic cleaner is mainly used in medical industry, computers, microelectronics computers. It uses high frequency to clean the sub micron dirt without damaging the equipment  
Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-506 HIGH FREQUENCY DESKTOP ULTRASONIC CLEANER



Digital display and control ultrasonic power adjustable range from 0 to 100%  
Special high frequency ultrasonic transducer (53/68/80/100/132 KHz)  
No stud- welding technology, Cleaning effect is better  
MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful  
MHC heating system, maximum temperature is up to 80°C  
All stainless steel material inner tank is made by punching SUS304 1.0mm  
2500V high voltage test to ensure the product safety, strict production, aging production quality control

## SPECIFICATIONS

Model	BULC-506
Capacity	15 L
Working frequency	53/68/80/100/132 KHz
Ultrasonic Power	360 W
Heating power	400 W/220V 400 W/110V
Timing time	1- 99 mins
Temperature range	0- 80°C
Power regulation	0-100 %
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	330x300x150 mm
Packing size	330x300x150 mm
Gross weight	11.5 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)