



SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-704

SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-704

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-704 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



- Standard basket and lid
- Time and temperature real time display
- Fully microprocessor controller and programmable
- Digital display, easy operation for time and temperature
- Low level water protection
- Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

SPECIFICATIONS

Model	BULC-704
Capacity	6 L
Ultrasonic frequency	40 KHz
Time range/ Heating temp	1-30 mins / 20- 80°C
Ultrasonic power	180 W
Heating power	400 W
Tank size	300x150x150 mm
Packing Size	430x280x380 mm
Gross Weight	6 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com