



# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-709

## SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-709

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

### BULC-709 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

## SPECIFICATIONS

Model	BULC-709
Capacity	22 L
Ultrasonic frequency	40 KHz
Time range/ Heating temp	1-30 mins / 20- 80°C
Ultrasonic power	480 W
Heating power	800 W
Tank size	500x300x150 mm
Packing Size	630x430x390 mm
Gross Weight	13 Kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)