



SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-712

SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-712

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-712 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



- Standard basket and lid
- Time and temperature real time display
- Fully microprocessor controller and programmable
- Digital display, easy operation for time and temperature
- Low level water protection
- Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

SPECIFICATIONS

| | |
|--------------------------|----------------------|
| Model | BULC-712 |
| Capacity | 2 L |
| Ultrasonic frequency | 80 KHz |
| Time range/ Heating temp | 1-480 mins / RT-80°C |
| Ultrasonic power | 50 W |
| Heating power | 100 W |
| Tank size | 150x140x100 mm |
| Packing Size | 260x260x280 mm |
| Gross Weight | 3.5 Kg |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com