



SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-712

SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-712

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-712 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

SPECIFICATIONS

Model	BULC-712
Capacity	2 L
Ultrasonic frequency	80 KHz
Time range/ Heating temp	1-480 mins / RT-80°C
Ultrasonic power	50 W
Heating power	100 W
Tank size	150x140x100 mm
Packing Size	260x260x280 mm
Gross Weight	3.5 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com