



# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-713

## SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-713

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

### BULC-713 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



- Standard basket and lid
- Time and temperature real time display
- Fully microprocessor controller and programmable
- Digital display, easy operation for time and temperature
- Low level water protection
- Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

## SPECIFICATIONS

Model	BULC-713
Capacity	4 L
Ultrasonic frequency	80 KHz
Time range/ Heating temp	1-480 mins / RT-80°C
Ultrasonic power	100 W
Heating power	400 W
Tank size	300x150x100 mm
Packing Size	445x250x400 mm
Gross Weight	8.5 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)