



## DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER BULC-803

## DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER BULC-803

Double adjustable frequencies ultrasonic cleaner uses frequency of 45 KHz /80 KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-803 DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER



- Over- Voltage and over current protection
- Standard basket and lid
- Low level water protection
- Fully microprocessor controller and programmable
- Time and temperature real time display
- Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

## SPECIFICATIONS

| Model                           | BULC-803            |
|---------------------------------|---------------------|
| Capacity                        | 27 L                |
| Adjustable ultrasonic frequency | 45 KHz /80 KHz      |
| Time range/ Heating temp        | 1-480 mins /0- 80°C |
| Ultrasonic power                | 500 W               |
| Heating power                   | 800 W               |
| Tank size                       | 500x300x180 mm      |
| Packing Size                    | 1100x690x800 mm     |
| Gross Weight                    | 101 Kg              |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)