



ULTRASONIC CLEANER BULC-905

ULTRASONIC CLEANER BULC-905

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-905 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion
LED display

SPECIFICATIONS

Model	BULC-905
Capacity	30 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	600 W
Heating power	500 W
Tank size	500x300x200 mm
External Size	530x330x350 mm
Packing Size	640x400x400 mm
Gross weight	15 Kg
Power Supply	AC 220V, 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com