



ULTRASONIC CLEANER BULC-907

ULTRASONIC CLEANER BULC-907

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-907 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion
LED display

SPECIFICATIONS

Model	BULC-907
Capacity	10 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-60 min/RT-80°C
Ultrasonic power	240 W
Heating power	500 W
Tank size	300x240x150 mm
External Size	330x270x310 mm
Packing Size	430x340x350 mm
Gross weight	8 Kg
Power Supply	AC 220V, 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com