



## ULTRASONIC CLEANER BULC-908

## ULTRASONIC CLEANER BULC-908

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-908 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe  
The shell and inner cavity are made of 304 stainless steel, which is anti corrosion  
LED display

## SPECIFICATIONS

Model	BULC-908
Capacity	15 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-60 min/RT-80°C
Ultrasonic power	360 W
Heating power	500 W
Tank size	330x300x150 mm
External Size	360x330x310 mm
Packing Size	460x400x400 mm
Gross weight	10 Kg
Power Supply	AC 220V, 50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)