



ULTRASONIC CLEANER BULC-909

ULTRASONIC CLEANER BULC-909

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-909 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe
The shell and inner cavity are made of 304 stainless steel, which is anti corrosion
LED display

SPECIFICATIONS

Model	BULC-909
Capacity	22.5 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-60 min/RT-80°C
Ultrasonic power	480 W
Heating power	500 W
Tank size	500x300x150 mm
External Size	530x330x310 mm
Packing Size	640x400x400 mm
Gross weight	13 Kg
Power Supply	AC 220V, 50Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com