



ULTRASONIC CLEANER BULC-911

ULTRASONIC CLEANER BULC-911

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-911 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe
The shell and inner cavity are made of 304 stainless steel, which is anti corrosion
LED display

SPECIFICATIONS

Model	BULC-911
Capacity	6.5 L
Ultrasonic frequency	Double frequency 28/40 Khz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	180 W
Heating power	500 W
Tank size	300x150x150 mm
External Size	330x180x310 mm



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com