



## ULTRASONIC CLEANER BULC-915

## ULTRASONIC CLEANER BULC-915

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-915 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion  
LED display

## SPECIFICATIONS

|                          |                            |
|--------------------------|----------------------------|
| Model                    | BULC-915                   |
| Capacity                 | 30 L                       |
| Ultrasonic frequency     | Double frequency 28/40 Khz |
| Time range/ Heating temp | 1-99 min/RT-80°C           |
| Ultrasonic power         | 600 W                      |
| Heating power            | 500 W                      |
| Tank size                | 500x300x200 mm             |
| External Size            | 530x330x350 mm             |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)