



## ULTRASONIC CLEANER BULC-916

## ULTRASONIC CLEANER BULC-916

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-916 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion  
LED display

## SPECIFICATIONS

Model	BULC-916
Capacity	6.5 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	180 W
Heating power	500 W
Tank size	300x150x150 mm
External Size	300x180x310 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)