



## MINI MIXER BVMX-302

## MINI MIXER BVMX-302

This sturdy, reliable mixers can adjust speed to meet different experimental requests. Multiple mixing modes, Touch and continuous operations. Aluminium cast base avoids unnecessary "walking". A variety of optional adapters are available which process 70 samples at the same time to make the experiment more convenient and efficient. Wide speed range:0-3000rpm. Built-in short mix, Microprocessor control, LED display the speed and time.

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

## BVMX-302 MINI MIXER



Mini and compact appearance, powerful function, sucker type under-chassis, super shockproof, ideal for high speed operation

Multiple mixing modes, touch-operation, continuous operation and timing set

Speed adjustable, and wide speed range, Stable low speed, powerful high speed

Shaking blocks convenient for replacement, stable and reliable, sturdy and durable eccentric bearing designed

Various tube holders are for optional (Suitable for Eppendorf Tubes)

## SPECIFICATIONS

Model	BVMX-302
Speed Range	0 ~ 2500 rpm
Shaking Orbit	3 mm
Shaking Mode	Circle
Operation	knob + scale
Operation Mode	continuous or press running
Voltage	AC 220V / AC 110V, 50/60 Hz
Power	60 W
Fuse	250 V, 1 A, $\Phi$ 5x20
Dimension	W.132 x D.156 x H.150 mm
Net Weight	2.5 kgs

## ACCESSORIES

Accessory Code	Name	Description
3100806020	Standard head	Suit for tubes and containers with dia.<30mm
3100806021	Standard platform	Standard platform to hold foam tube holder

## OPTIONAL ACCESSORIES

Accessory Code	Name
3100806006	Convex pad
3100806007	Foam tube holder

Accessory Code	Name
3100806013	Tube holder platform
3100806014	Tube holder platform

3100806008	Foam tube holder
3100806009	Foam tube holder
3100806010	Foam tube holder
3100806011	Foam tube holder
3100806012	Tube holder platform

3100806015	Plate holder platform
3100806016	Universal platform for tube stands and all standard plates
3100806017	Universal platform for tube stands and all standard plates
3100806018	Universal platform for tube stands and all standard plates
3100806019	Customized



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
 Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)