



COD ANALYZER

COD ANALYZER BFU1T1

COD ANALYZER



1. 7 inch high-resolution touch screen with friendly GUI.
2. The unique high-precision filter system makes the instrument more stable and accurate.
3. No blank calibration, automatic measurement, direct reading of concentration.
4. With built-in standard curves. Users can set up curves according to different measurement requirements.
5. Support zero point and up to 7-point calibration. Support calibration restoring function to restore factory calibration data.
6. With a large-capacity storage function, can store 5000 sets of test data.
7. Built-in thermal printers can print current data and stored data.

SPECIFICATIONS

Model	BFU1T1
Measurement Method	Fast digestion spectrophotometry
Measuring Range	10-10000mg/L
Accuracy	≤5%
Resolution	0.01
Light Source	Imported Silicon Photodiode
Display	3-inch color screen
Output	Micro USB
Power	Power adapter/Rechargeable lithium battery
Waterproof Level	IP54
Ambient Temperature	(5-40) °C
Ambient Humidity	Relative humidity ≤ 85% (No condensation)
Dimensions/N.W.	100x210x70 mm/1kg
Shipping size/G.W.	355x265x150mm/2kg
Standard Accessories	Sample bottle x 2 pcs, Power adapter(Charger), Manual
Alt Name	COD Analyzer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com