



PORTABLE STERILIZER

PORTABLE STERILIZER BER1CK1 BER1CK2

PORTABLE PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Electric or LPG heated.
3. Double scale indication pressure gauge.
4. Easy to operate, safe and reliable.

SPECIFICATIONS

Model	BER1CK1	BER1CK2
Capacity	18L	24L
Working Pressure	0.14-0.16Mpa	
Working Temperature	126°C	
Max.Safety Pressure	0.165Mpa	
Power	AC220V.50Hz/2Kw	
Packing Size(mm)	410*410*430	410*410*550
G.W/N.W	16/14KG	17/15KG
Alt Name	Portable Pressure Steam Sterilizer	

PORTABLE STERILIZER BER1CJ1 BER1CJ2

PORTABLE PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Double scale indication pressure gauge.
3. Electric heated.
4. Easy to operate, safe and reliable.

SPECIFICATIONS

Model	BER1CJ1	BER1CJ2
Capacity	18LΦ280*260mm	24LΦ280*390mm
Working Pressure	0.14-0.16Mpa	
Working Temperature	126°C	

Max.Safety Pressure	0.165Mpa	
Heat Average	$\leq \pm 0.1^{\circ}\text{C}$	
Power	2KW/AC200V.50Hz	
Packing Size	410*410*520mm	410*410*630mm
G.W/N.W	18kg/16kg	19kg/17kg
Alt Name	Portable Pressure Steam Sterilizer	

PORTABLE STERILIZER BER1CL1 BER1CL2

PORTABLE PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Hand wheel type of quick-open door structure.
3. Indicator light indicates working state.
4. Double scale indication pressure gauge.
5. Time scope: 0~60 min.
6. Automatically shut off with beep reminding after sterilization.
7. heated.
8. Easy to operate, safe and reliable.

SPECIFICATIONS

Model	BER1CL1	BER1CL2
Capacity	18L	24L
Working Pressure	0.14-0.16Mpa	
Working Temperature	126°C	
Temp. range	126°C	
Max.Safety Pressure	0.165Mpa	
Timer range	0-60min	
Heat average	$\leq \pm 1^{\circ}\text{C}$	
Power	AC220V.50Hz/2Kw	
Packing Size(mm)	410*410*430	410*410*550
G.W/N.W	16/14KG	17/15KG
Alt Name	Portable Pressure Steam Sterilizer	

PORTABLE STERILIZER BER1CM1 TO BER1CM4

VERTICAL PRESSURE STEAM STERILIZER



1. Fully stainless steel structure.
2. Hand wheel type of quick-open door structure.
3. Door safety lock system.
4. Digital display of working status, touch of key.
5. Self-inflating type seal.
6. Over temperature & over pressure auto-protection.
7. Safe protection of water lacking.
6. Automatically shut off with beep reminding after sterilization.
7. Auto discharge the cool air, and steam discharging automatically after sterilization.
8. Used for sterilization of medical instruments, medical cotton products.
9. Fully stainless steel SUS304 / AISI 304 - 3mm.

SPECIFICATIONS

Model	BER1CM1	BER1CM2	BER1CM3	BER1CM4
Capacity	35L(Φ318*450)mm	50L(Φ340*550)m m	75L(Φ400*600)m m	100L(Φ440*650)m m
Working pressure	0.22Mpa			
Working temperature	134℃			
Max.working pressure	0.23Mpa			
Heat average	≤± 1℃	≤ ± 1℃		≤± 1℃
Timer	0-99min or 0-99hour59 min			
Temp.range	105℃~134℃			
Power	2.5KW/AC220V.50Hz	3KW/AC220V.50Hz	4.5KW/AC220V.50Hz	
Exterior size	450*450*1010mm	510*470*1130mm	560*560*1120mm	540*560*1250mm
Packing size	570*550*1150mm	590*590*1280mm	650*630*1280mm	680*630*1370mm
G.W/N.W	72kg/56kg	88kg/68kg	100kg/80kg	110kg/85kg
Alt Name	Vertical Pressure Steam Sterilizer			



FEATURES

1. Fully stainless steel structure.
2. Hand wheel type of quick-open door structure.
3. Door safety lock system.
4. Digital display of working status, touch of key.
5. Self-inflating type seal.
6. Over temperature & over pressure auto-protection.
7. Safe protection of water lacking.
6. Automatically shut off with beep reminding after sterilization.

7. Auto discharge the cool air, and steam discharging automatically after sterilization.
8. Used for sterilization of medical instruments, medical cotton products.
9. Fully stainless steel SUS304 / AISI 304 - 3mm.
10. Possible to install drying system according request.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com