



HEMATOCRIT CENTRIFUGE

HEMATOCRIT CENTRIFUGE BKA-2501

MICRO HEMATOCRIT AND MICRO TUBES CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,Program(Memory)information, Acceleration& 3.Deceleration information and Error code information.
4. ABS Centrifuge body and in light weight 8KGS.
5. Safety Electrical lid interlock,machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuous according the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mentioned when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the scre

SPECIFICATIONS

Model	BKA-2501
Max. Speed	16000rpm
Max. RCF	17800xg
Max. Capacity	24 Placers
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	460*350*265mm
Net weight(without rotor)	8KGS
Power supply	AC110~240V,50/60HZ,3A,0.2KW
Alt Name	Micro Hematocrit and Micro Tubes Centrifuge

ACCESSORIES FOR PURCHASE

No	name	Capacity (ml)	Max. Speed (rpm)	Max. RCF (xg)
1	NO.1 HEMATOCRIT ROTOR	24 placers	12000rpm	14800xg
2	NO.2 MICRO TUBE ROTOR	12x1.5/2.2ml	16000rpm	17800xg

FEATURES

1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,Program(Memory)information, Acceleration& 3.Deceleration information and Error code information.
4. ABS Centrifuge body and in light weight 8KGS.
5. Safety Electrical lid interlock,machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuous according the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mentioned when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the screen.
11. 24 placers hematocrit rotor and micro rotor for choice.
12. With CE and ISO9001:2008,ISO13485 certificate

HEMATOCRIT CENTRIFUGE BKA-1201

HEMATOCRIT CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed, timer, RCF, rotors information, Program (Memory) information, Acceleration & Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant to detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock, door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.

SPECIFICATIONS

Model	BKA-1201
Max. Speed	16000rpm
Max. RCF	14800xg
Max. Capacity	24PCS
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	330*435*340mm
Net weight(without rotor)	28KGS
Power supply	AC220-240V, 50/60HZ, 5A, 0.5KW
Alt Name	Hematocrit Centrifuge

ACCESSORIES FOR PURCHASE

No	name	Capacity (mL)	Max Speed (rpm)	Max RCF (xg)
1	capillary rotor	24 pcs capillary rotor	12000rpm	14800xg

FEATURES

1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed, timer, RCF, rotors information, Program (Memory) information, Acceleration & Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant to detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock, door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.
9. 10 kinds of Program (Memory) and 10 kinds of Acceleration & Deceleration for choice. If ask for more Programs please specially mention when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the screen.
11. Wide choice of Rotors and adapters.
12. With CE and ISO9001:2008, ISO13485 certificate.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com