



CYTO CENTRIFUGE

CYTO CENTRIFUGE BKA-2101

CYTO(TCT) CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.

SPECIFICATIONS

Model	BKA-2101
Max. Speed	2200rpm
Max. RCF	860xg
Max. Capacity	12 samples
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	540*430*265mm
Net weight(without rotor)	28KGS
Power supply	AC110~240V,50/60HZ,1A,0.04KW
Alt Name	Cyto(TCT) Centrifuge

ACCESSORIES FOR PURCHASE

No	name	Capacity (ml)	Max. Speed (rpm)	Max. RCF (xg)
1	NO.1	double layer total 8 samples	2200rpm	395xg
2	NO.2	double layer total 6 samples	2200rpm	735xg
3	NO.3	12 samples	2200rpm	860xg

FEATURES

1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mention when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the screen.
11. Wide choice of Rotors and adapters.
12. With CE and ISO9001:2008,ISO13485 certificate.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com