



GEL CARD CENTRIFUGE

GEL CARD CENTRIFUGE BCBLS-106

BLOOD GROUP SEROLOGY CENTRIFUGE




SPECIFICATIONS

Model	BCBLS-106
Max. Speed	3840rpm
Max. RCF	1790xg
Max. Capacity	24 gel cards
Drive System	Direct, induction motor drive
Noise Level	≤ 50 dB(A)
Power Supply	AC 200-240 V, 50/60 Hz or AC 100-127 V, 50/60 Hz
Max. Power Consumption	100 W
Dimension (W x D x H)	14.8 x 14.2 x 11.8 in / 37.5 x 36 x 30 cm
Weight w/o Rotor	46 lb / 21 kg
Alt Name	Gel Card Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Rotor Capacity
1	No. 1 Swing-out Rotor	12 gel cards
2	No. 2 Swing-out Rotor	24 gel cards



GEL CARD CENTRIFUGE BIFG-502 BIFG-503

GEL CARD CENTRIFUGE

The gel card centrifuge holds automatic balance control and the user can set speed and centrifugation time. It can be suitable for many types of gel cards from different manufacturers.



Digital interface visor control and visualization of centrifuge speed and time remaining for processing.

Automatic balance control.

Pre-defined centrifugal parameters for optimum gel cards procession, acceleration, speed, braking and time.

Programmable audible alarms for end of centrifuge time periods.

SPECIFICATIONS

Model	BIFG-502	BIFG-503
Max. Speed	1800rpm	1500rpm
Max. RCF	272xg	294xg
Max. Capacity	12 gel cards	24 gel cards
Timer	1~20 min	
Noise	≤55 db(A)	
Power	AC 100-240 V 85 W	
Dimension (LxWxH) mm	440x340x170 mm	530x410x180 mm
Weight	7 kg	11.5 kg
Alt Name	Gel Card Centrifuge	



BIFG-502



BIFG-503

GEL CARD CENTRIFUGE BCBL-308

CENTRIFUGE FOR GEL CARDS



1. Brushless DC motor to make sure free maintenance no powder pollution, micro computer control and LCD display the time, speed and RCF to let operation easily.
2. There is specialized step-ladder centrifugation technology; different step-ladders could be designed by different reagents.
3. There are several kinds of regular centrifugation programs; other centrifugation programs also could be designed by customers' demand.
4. Electronic lock, over speed and imbalance protection, there is a self-check protection program, safe and reliable, simple operation, scientific and practicable.
5. With CE & ISO9001 & ISO13485 Approval.

SPECIFICATIONS

Model	BCBL-308
Max. Speed	3000rpm
Max. RCF	1190xg
Max. Capacity	24 cards/12 cards
Noise	≤55dBA
Timer	0~99min59s
Power supply	AC 220 50HZ 2A
Dimension (12 cards)	460*320*230mm
Net. Weight (12 cards)	17Kgs
Dimension (24 cards)	500*360*250mm
Net. Weight (24 cards)	19Kgs
Speed accuracy	±20rpm
Package	carton box
Alt Name	Centrifuge for Gel Cards

APPLICATIONS

It is widely used in the test of Blood Group Identification,Cross-Matching of Blood,Irregular Antibody Screening etc.

GEL CARD CENTRIFUGE BKA-2301

BLOOD ID CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.

SPECIFICATIONS

Model	BKA-2301
Max. Speed	2200rpm
Max. RCF	2200xg
Max. Capacity	12 ID
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)	460*350*265mm
Net weight(without rotor)	23KGS
Power supply	AC110~240V,50/60HZ,1A,0.04KW
Alt Name	Blood ID Centrifuge

FEATURES

1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially mention when place the order.
10. Automatic troubleshooting with indication of error code and alarm on the screen.
11. Wide choice of Rotors and adapters.
12. With CE and ISO9001:2008,ISO13485 certificate.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com