



BULLET ICE MAKER

BULLET ICE MAKER BFC1C1 TO BFC1C5

VERTICAL ICE MAKER



Luxurious appearance and compact design;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Cylindrical bullet ice with elegant appearance.
Ice thickness can be adjusted, with a wider scope of application.
Bullet Ice

SPECIFICATIONS

Model	BFC1C1	BFC1C2
Old Model	BIBU-101	BIBU-102
Ice making capacity	20kg/24h	25kg/24h
Power/voltage	220w	240w
Ice storage capacity	6kg	7kg
Ice shape	Bullet shape	
N.W./kgs	28kg	30kg
G.W./kgs	31kg	32kg
Product dimension	330x470x605mm	398x542x682mm
Packing Dimension	400x550x680	480x630x730mm
20GP	168	135
40GP	348	285
40HQ	464	285
Alt Name	Vertical Ice Maker	

Model	BFC1C3	BFC1C4	BFC1C5
Old Model	BIBU-104	BIBU-109	BIFL-212
Ice making capacity	50kg/24h	80kg/24h	120kg/24h
Power/voltage	380w	580w	760w
Ice storage capacity	14kg	25kg	50kg
Ice shape	Bullet shape		
N.W./kgs	35kg	51kg	70kg
G.W./kgs	38kg	56kg	75kg
Product dimension	400x510x825mm	500x612x900mm	680x610x1050mm
Packing Dimension	480x590x890mm	580x700x970mm	750x680x1100mm
20GP	98	64	48
40GP	196	136	96
40HQ	290	136	96
Alt Name	Vertical Ice Maker		



BFC1C1



BFC1C2



BFC1C3



BFC1C4



BFC1C5

FEATURES

Luxurious appearance and compact design;

Fully automatically computer controlled, convenient for use;

Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;

Highly efficient fluoride-free compressor, energy efficient and environment-friendly;

Cylindrical bullet ice with elegant appearance.

Ice thickness can be adjusted, with a wider scope of application.



Bullet Ice



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com