



## DELIVERY TABLE

## DELIVERY TABLE BGY1H1

### ELECTRIC GYNECOLOGICAL TABLE

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Whole stainless steel hidden auxiliary table with flexible extension and retraction  
Equipped with hidden sewage basin  
With foot control, can adjust the movements easily  
Imported LINAK gear motor, without noise, have a long service life

## SPECIFICATIONS

Model	BGY1H1
Old Model	BHBD-211
Table Length	1450+50 mm
Table width	520+30 mm
Height Adjustment	(400 900)+50 mm
Trendelenburg/rev-trend	≥15°/≥25°
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin
Power Supply	220V 50Hz
Alt Name	Electric gynecological table

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation

## DELIVERY TABLE BGX1D1

### MANUAL MULTI-FUNCTION GYNAECOLOGICAL TABLE

The Delivery table is a flexible design suitable for conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Material: Steel Powder Coated S.S, Color: Blue

## SPECIFICATIONS

Model	BGX1D1
Old Model	BHBD-206
Function	Multi
Back-rest lifting angle	≥20°
Back section up adjusted	≥60°
Seat board upwards	≥8°
Seat board downwards	≥10°
Load Bearing	250 kg
Height adjustment	≤650 mm
Tabletop Dimension	1240Lx600W mm
Upright Dimension	500Lx600W mm
Dimension	1750Lx600Wx650~900 mm
Alt Name	Manual Multi-Function Gynaecological Table

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation.

## DELIVERY TABLE BGY1G1

### ELECTRIC GYNECOLOGICAL TABLE

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Adopt high density memory sponge and is once formed, with no seam  
Hydraulic lifting by foot pedal, other actions are controlled by pneumatic  
Imported Y shaped sealing ring and durable  
Stainless steel concealed auxiliary table  
The base is made up of high strength ABS  
Equipped with foot brake device which is easy for moving the table

## SPECIFICATIONS

Model	BGY1G1
Old Model	BHBD-210
Table Length	1300+50 mm
Table width	550+30 mm
Height Adjustment	(700 950)+50 mm
Trendelenburg/rev-trend	≥3°/≥10°
Dimension of auxiliary platform	(350x450)+20 mm
Backboard downward	≥70°
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin
Alt Name	Electric gynecological table

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation

## DELIVERY TABLE BGY1C1

### ELECTRIC GYNECOLOGICAL TABLE

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Various actions are controlled by Imported brand motor  
Hand-held controller and foot controller, mutually independent and is easily for operation  
Leg support frame can fold up and down electric  
Unique profile design which is beautiful and elegant  
Whole stainless steel hidden auxiliary table with flexible extension and retraction  
Equipped with hidden sewage basin  
The table adopts high density memory sponge and is once foam molding, with no seam  
Features with one-key reset function  
Super mute foot wheel and equipped with brake device

## SPECIFICATIONS

Model	BGY1C1
Old Model	BHBD-204
Back board length	900 mm
Seat board length	400 mm
Table width	600 mm
Trendelenburg/rev-trend	$\geq 10^\circ/\geq 12^\circ$
Backboard Adjustment	$\geq 45^\circ$
Height adjustment (from back top to ground)	1500-1750 mm
Standard Accessory	A pair of arm rods, a pair of leg support, a stainless steel sewage basin, auxiliary table (already installed), a foot controller
Alt Name	Electric Gynecological table

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation

## DELIVERY TABLE BGY1E1

### DELIVERY TABLE 4002

The Delivery table is a flexible design suitable for conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Imported LINAK motor various actions are controlled by foot controller  
The base seat and upright column covering are made of high quality 304 stainless steel which is anti corrosive and can be easily clean  
External folding auxiliary table is for easy dismantling  
Equipped with stainless steel sewage basin and other accessories  
Foot treadle brake device which is easy for moving obstetric table

## SPECIFICATIONS

Model	BGY1E1
Old Model	BHBD-208
Table Length	2000+50 mm
Table width	600+10 mm
Height Adjustment	(650-900)+50 mm
Trendelenburg/rev-trend	$\geq 10^\circ / \geq 12^\circ$
Backboard upward	$\geq 50^\circ$
Leg board outward	$\geq 90^\circ$
Working arealength/width	600x600 mm
Standard accessory	A pair of arm support, a pair of leg support, a pair of foot pedal, a pair of foot arm rods, a stainless steel sewage basin, a foot controller
Alt Name	Delivery table 4002

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation

## DELIVERY TABLE BGY1F1

### ELECTRIC GYNECOLOGICAL TABLE

The Delivery table is a flexible design for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.



Various actions are controlled by Imported LINAK motor  
Hand-held remote control and pedal remote control, independent of each other, easy to operate it  
With concealed stainless steel sewage waste basin  
The mattress is wrapped with imported artificial leather. High density memory sponge is formed in one piece, with no seam  
Equipped with brake device which is easy for moving the table  
The base is made up of high strength ABS, which is corrosion resistant and easy to clean  
Leg plate can be lifted electrically  
Concealed stainless steel auxiliary table  
With function of auto restoration

## SPECIFICATIONS

Model	BGY1F1
Old Model	BHBD-209
Table Length	1270+50 mm
Table width	550+30 mm
Height Adjustment	(700-1000)+50 mm
Trendelenburg/rev-trend	≥5°/≥12°
Backboard downward	≥40°
leg board reflex	≥30°
Working area	(350x480)+20 mm
Standard accessory	A pair of arc shaped handle, a pair of leg support, a stainless steel sewage basin, a power cord, a hand held remote controller, a foot controller
Main Voltage AC	220V+22V, 50Hz+1Hz
Alt Name	Electric gynecological table

## APPLICATIONS

Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)