



PIPETTE PUMP

PIPETTE PUMP BFS1BH1 BFS1BH2 BFS1BH3

PIPETTE PUMP

We are engaged in offering a wide range of Pipette Pumps that provide safe and accurate pipetting and dispensing of liquids. It can be rotated to dispense and draw up the liquid. Available in a variety of sizes and shapes.



Economical
Highly gas tightness
Environment- friendly materials
Lightweight and dexterous for operation
Durable
Recalibration Mechanism

SPECIFICATIONS

Model	BFS1BH1	BFS1BH2	BFS1BH3
Volume	2ml	10ml	25ml
Colour	Blue + White	Green + White	Red + White
Unit	1000		
Alt Name	Pipette Pump		



BFS1BH1



BFS1BH2



BFS1BH3

APPLICATIONS

Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com