


## SEROLOGICAL PIPETTE

# SEROLOGICAL PIPETTE BFS1BN1 TO BFS1BN14

## SEROLOGICAL PIPETTE

Serological pipettes are designed to be used with electromotion pipette controllers and manual pipette pumps to transfer liquid in various volumes for laboratory use. With high transparency polystyrene and clear graduation, the serological pipettes are easy to read in accurate volume by users. Graduation is printed on the pipette surface to be calibrated on accuracy  $\pm 2\%$ . To meet global market needs, pipette tips are crafted to two forms, welding and stretching.

Two types of packaging are available for pipettes: bulk packaging and individually packaging in plastic/paper wrapper. All serological pipettes are packaged with sterilization in bags, RNase-free, DNase-free and Non-pyrogenic.



Easy to use  
Medical grade original Polystyrene  
RNase- free, DNase- free and Non-pyrogenic  
Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml  
Bidirectional graduation improve accuracy  
Seven types of color code on tips make pipette  
Sterilized by Gamma irradiation

## SPECIFICATIONS

Model	BFS1BN1	BFS1BN2	BFS1BN3	BFS1BN4	BFS1BN5
Old Model	BPIC-601	BPIC-602	BPIC-603	BPIC-604	BPIC-605
Volume	1 ml		2 ml		5 ml
Precision	2%				
Negative Graduation (ml)	-0.4				
Graduation	1 / 100				10-Jan
Color	Yellow				Blue
Sterile	Gamma				
Package	Individual	Bulk pack			Individual
Pcs/bag	100 pcs/bag		50 pcs/bag	25 pcs/bag	100 pcs/bag
Pcs/ctn	4000 pcs/carton	3200 pcs/carton	1800 pcs/carton	1200 pcs/carton	800 pcs/carton
Unit	1				
Alt Name	Standard Serological Pipette				

Model	BFS1BN6	BFS1BN7	BFS1BN8	BFS1BN9	BFS1BN10
Old Model	BPIC-606	BPIC-607	BPIC-608	BPIC-609	BPIC-610
Volume	5 ml	10ml		25ml	
Precision	2%				
Negative Graduation (ml)	-0.4	-3			
Graduation	10-Jan			05-Jan	
Color	Blue	Red		Purple	
Sterile	Gamma				
Package	Bulk pack	Individual	Bulk pack	Individual	Bulk pack
Pcs/bag	25 pcs/bag				
Pcs/ctn	600 pcs/carton	4000 pcs/carton	3200 pcs/carton	3000 pcs/carton	2100 pcs/carton
Unit	1				
Alt Name	Standard Serological Pipette				

Model	BFS1BN11	BFS1BN12	BFS1BN13	BFS1BN14
Old Model	BPIC-611	BPIC-612	BPIC-613	BPIC-614
Volume	50 ml		100 ml	
Precision	2%			
Negative Graduation (ml)	-3		-	

Graduation	1			
Color	Black			
Sterile	Gamma			
Package	Individual	Bulk pack	Individual	Bulk pack
Pcs/bag	25 pcs/bag		1 pcs/bag	
Pcs/ctn	1400 pcs/carton	1000 pcs/carton	300 pcs/carton	
Unit	1			
Alt Name	Standard Serological Pipette			

## APPLICATIONS

Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)