



ETHYLENE OXIDE STERILIZER

ETHYLENE OXIDE STERILIZER BKM-401 TO BKM-407

ORDINARY ETHYLENE OXIDE STERILIZER



The series of epoxyethane sterilizer is a specialized sterilization equipment that uses a mixture of pure epoxyethane and epoxyethane gas to fumigate items enclosed in a sterilization room at a certain temperature, pressure, humidity, and duration. Ethylene oxide gas is a very active chemical disinfectant, which not only achieves broad-spectrum sterilization but also reliable sterilization.

SPECIFICATIONS

Model	BKM-401	BKM-402	BKM-403	BKM-404	BKM-405	BKM-406	BKM-407
Inner Gallbladder Capacity	72L	120L	170L	195L	271L	304L	508L
Disinfection Mode	Low temperature						
Source	AC220V 50Hz						
Inner Liner Size (L x W x H) mm	60 x 30 x 40	60 x 40 x 50	67 x 43 x 59	70 x 45 x 62	71.5 x 55 x 69	80 x 55 x 69	96 x 67 x 79
Outer Box Size (L x W x H) mm	78 x 51 x 113	78 x 61 x 123	85 x 64 x 132	88 x 66 x 135	89.5 x 76 x 142	98 x 76 x 142	118 x 88 x 160
Alt Name	Ordinary Ethylene Oxide Sterilizer						

ETHYLENE OXIDE STERILIZER BKM-101 TO BKM-108

FULLY AUTOMATIC ETHYLENE OXIDE STERILIZER



The fully automatic ethylene oxide sterilizer is controlled by a microcomputer color LCD and is divided into a computer control system, a vacuum sterilization room system, a constant temperature system, and a residual gas treatment system. When in use, you only need to press the run button twice to achieve automatic heating, vacuuming, dosing, humidification, vaporization, sterilization, residue removal, and automatic analysis of the ethylene oxide residue. In order to make our product more perfect, we have added single tank pure ethylene oxide dosing and mixed ethylene oxide dosing functions on this machine for more convenience and safety, achieving 100% sterilization. This product has unique advantages in sterilization and disinfection of moisture and heat sensitive items.

SPECIFICATIONS

Model	BKM-101	BKM-102	BKM-103	BKM-104	BKM-105	BKM-106	BKM-107	BKM-108
Inner Gallbladder Capacity	61L	111L	130L	163L	182L	203L	301L	400L
Disinfection Mode	Low temperature							

Source	AC220V 50Hz							
Inner Liner Size (L x W x H, mm)	620 x 350 x 280	610 x 480 x 380	710 x 480 x 380	770 x 460 x 460	860 x 460 x 460	960 x 460 x 460	960 x 560 x 560	1010 x 660 x 600
Outer Box Size (L x W x H, mm)	770 x 660 x 1400	760 x 790 x 1600	860 x 790 x 1600	920 x 700 x 1760	1010 x 770 x 1760	1110 x 770 x 1760	1110 x 870 x 1960	1160 x 1130 x 1950
Alt Name	Fully Automatic Ethylene Oxide Sterilizer							



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com