



SAMPLE TUBE SCANNER

SAMPLE TUBE SCANNER BFS1BS1 BFS1BS2

SAMPLE TUBE SCANNER

Sample tube scanner is specifically designed to operate on top of the sample tubes without applying any weight, avoiding any damage.



It can be finished within 5 seconds for one of 96 holes rack or 384 holes rack. Advanced decoder technology can be identified as stable and secure even if the code is fuzzy.

Optional item for 1-D code rack scanning function.

It can scan the difference of height of the sample tube in the same rack .

Compact, ergonomic design assures efficiency of workflow.

Detection of rack type automatically for setting optimum scanning.

SPECIFICATIONS

Model	BFS1BS1	BFS1BS2
Rack type available	13,24,48,96,384 holes rack	
1-D Code available	2/5 Industrial / Interleaved, Code 39, Code 128, Pharmacode, Codabar, EAN 13	
2-D Code available	Datamatrix ECC 200, PDF417, QR Code	
Camera	10 mega-pixel camera CMOS	
Alt Name	Sample tube scanner	

APPLICATIONS

Laboratory, Research, Medical



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com