





ALL STEEL BENCH


ALL STEEL BENCH BKP-201 BKP-202

STEEL WOOD BENCH




SPECIFICATIONS

Model	BKP-201	BKP-202
Size(W*D*H)	1000*750*850mm	1500*750*850mm
Countertop	12.7mm thick Solid Chemical-Resistant Phenolic Resin Board	
Cabinets	1 Cabinets	2 Cabinets
Structure	C-frame	
Drawer Slides	16-inch DTC brand dedicated silent slide rails	
Hinges	DTC brand hinges, 115-degree opening angle	
Handles	Aluminum alloy U-shaped handles or PP embedded handles	
Leveling Feet	Adjustable nylon base, acid and alkali resistant	
Alt Name	Steel Wood Bench	




BKP-201



BKP-202


ALL STEEL BENCH BKN-101 BKN-102 BKN-103

STEEL BENCH




SPECIFICATIONS


Model	BKN-101	BKN-102	BKN-103
Size (W*D*H)	1000*750*850mm	1500*750*850mm	
Structure	All-steel structure		
Countertop	12.7 mm physicochemical board countertop	12.7mm physicochemical board	
Hinges	DTC 105° hinges		
Drawer Slides		Silent drawer slides	
Sink			PP sink
Faucet			Three-way faucet
Design	As shown in images		
Alt Name	Steel Bench		



BKN-101




BKN-102



BKN-103


ALL STEEL BENCH BKP-101 BKP-102

ALL STEEL BENCH




SPECIFICATIONS

Model	BKP-101	BKP-102
Size(W*D*H)	1500*750*850mm	
Countertop	12.7mm thick Solid Chemical-Resistant Phenolic Resin Board	
Cabinet	2 cabinets	2 Cabinets
Water Sink		1 Water Sink
Water Tap		1 water tap
Alt Name	All Steel Bench	




BKP-101



BKP-102


ALL STEEL BENCH BKN-301 TO BKN-304

STEEL BENCH




SPECIFICATIONS


Model	BKN-301	BKN-302	BKN-303	BKN-304
Size (W*D*H)	1000*750*850mm	1500*750*850mm	1000*750*850mm	1500*750*850mm
Structure	All-steel structure			
Countertop	12.7 mm physicochemical board countertop			
Hinges	DTC 105° hinges			
Drawer Slides	Silent drawer slides			
Design	As shown in images			
Alt Name	Steel Bench			




BKN-301



BKN-302



BKN-303



BKN-304

ALL STEEL BENCH BKJ-101 BKJ-102

ALL STEEL LAB SIDE BENCH




SPECIFICATIONS


Model	BKJ-101	BKJ-102
Size (W*D*H)	1500*750*850mm	
Structure	All-steel, floor-standing structure	
Cabinet & Drawer Material	1.0 mm thick cold-rolled steel plate	
Surface Finish	Epoxy coated	
Table Top	12.7 mm thick black phenolic resin top	
Drawer Slides	Three-section silent drawer slides	
Hinges	Soft-close hinges	
Handles	All-steel straight pull handles	
Water Tap		1 water tap
Water Sink		1 water sink
Configuration	Two cabinets - one left and one right	1 basin and 1 single cabinet
Alt Name	All Steel Lab Side Bench	



BKJ-101



BKJ-101



BKJ-102



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com