





STAINLESS STEEL BENCH

STAINLESS STEEL BENCH BKN-201

STAINLESS STEEL BENCH



SPECIFICATIONS

Model	BKN-201
Size (W*D*H)	1500*750*850mm
Structure	Stainless steel structure
Countertop	12.7 mm physicochemical board countertop
Hinges	DTC 105° hinges
Drawer Slides	Silent drawer slides
Design	As shown in images
Alt Name	Stainless Steel Bench


STAINLESS STEEL BENCH BKJ-201

ALL STAINLESS STEEL LAB SIDE BENCH



SPECIFICATIONS

Model	BKJ-201
Size (W*D*H)	1500*750*850mm
Structure	All steel and floor-standing structure
Cabinet and drawer	1.0mm thick stainless steel plate
Table Top	12.7mm thick black phenolic resin top
Drawer Slides	Three-section silent drawer slides
Hinges	Soft-close hinges
Handle	All stainless steel straight pull handle
Configuration	Two cabinets - left and right cabinet
Alt Name	All stainless steel lab side bench



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com