



PP BENCH

PP BENCH BKJ-401 BKJ-402

ALL PP LAB SIDE BENCH




SPECIFICATIONS

Model	BKJ-401	BKJ-402
Size (W*D*H)	1500*750*850mm	
Structure	All PP, floor-standing structure	
Cabinet and Drawer	8 mm PP (polypropylene) sheet, welded using rods of the same color and material	
Table Top	12.7 mm thick black phenolic resin top	
Drawer Slides	PP drawer slides	
Hinges	PP hinges	
Handle	C-type or all-PP straight pull handle	
Water Tap	_____	1 water tap
Water Sink	_____	1 water sink
Configuration	Two cabinets - left and right cabinet	1 basin and 1 single cabinet
Alt Name	All PP lab side bench	



BKJ-401



BKJ-402

PP BENCH BKP-401 BKP-402

PP LAB BENCH




SPECIFICATIONS

Model	BKP-401	BKP-402
Size(W*D*H)	1500*750*850mm	
Countertop	8mm pure PP board (or optional 12.7mm phenolic resin board)	
Cabinets	2 Cabinets	
Cabinet Structure	Modular assembly structure with concealed internal connections, zero exposed metal screws	
Main Material	8mm thick high-grade pure PP (Polypropylene) panels, highly resistant to strong acids and alkalis	
Drawer Slides	PP self-lubricating sliding system (metal-free design)	
Hinges & Screws	PP material, 180-degree opening angle, highly corrosion-resistant	
Handles	Integrated U-shaped matte PP handles	
Alt Name	PP Lab Bench	




BKP-401



BKP-402


PP BENCH BKN-401 BKN-402

PP BENCH




SPECIFICATIONS

Model	BKN-401	BKN-402
Size (W*D*H)	1500*750*850mm	
Structure	PP structure	
Countertop	12.7mm physicochemical board countertop	
Sink	PP sink	
Faucet	Three-way faucet	
Design	As shown in images	
Alt Name	PP Bench	



BKN-401



BKN-402



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com