



DENTAL AUTOCLAVE BADT-101

DENTAL AUTOCLAVE BADT-101

AUTOCLAVE



- 1.Safety Valve
Prevent chamber from over pressure
- 2.Bacteria Filter
High quality,2µm precision
- 3.Steam Generator
Patented design, high efficiency and reliable
- 4.Air Switch
Shut down the machine in case of electricity leakage or short circuit
- 5.Double Fans
Bring faster cooling and more efficiency sterilization
- 6.Water Steam Separator
Patented design, prevent vacuum pump from oil damaging

SPECIFICATIONS

Model	BADT-101
Voltage	220V / 50Hz
Power	1750W
Capacity	8L
Chamber size	Φ17*32cm
Product size	53.5*36.5*32cm
Installation size	50cm
N.W.	35kgs
Printer	optional
Alt Name	Table-top autoclave



FEATURES



- 1.Safety Valve
Prevent chamber from over pressure



2. Bacteria Filter
High quality, 2µm precision



3. Steam Generator
Patented design, high efficiency and reliable



4. Air Switch
Shut down the machine in case of electricity leakage or short circuit



5. Double Fans
Bring faster cooling and more efficiency sterilization





6. Water Steam Separator

Patented design, prevent vacuum pump from oil damaging



7. Brass Piping

Withstand high temperature and pressure safe and reliable



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com