



## DENTAL AUTOCLAVE BADT-103

# DENTAL AUTOCLAVE BADT-103

AUTOCLAVE



- 1.Safety Valve  
Prevent chamber from over pressure
- 2.Bacteria Filter  
High quality,2µm precision
- 3.Steam Generator  
Patented design, high efficiency and reliable
- 4.Air Switch  
Shut down the machine in case of electricity leakage or short circuit
- 5.Double Fans  
Bring faster cooling and more efficiency sterilization
- 6.Water Steam Separator  
Patented design, prevent vacuum pump from oil damaging

## SPECIFICATIONS

Model	BADT-103
Voltage	220V / 50Hz
Power	1750W
Capacity	18L
Chamber size	Φ25*35cm
Product size	66*44.5*39cm
Installation size	60cm
N.W.	47kgs
Printer	optional
Alt Name	Table-top autoclave



## FEATURES



- 1.Safety Valve  
Prevent chamber from over pressure



2. Bacteria Filter  
High quality, 2µm precision



3. Steam Generator  
Patented design, high efficiency and reliable



4. Air Switch  
Shut down the machine in case of electricity leakage or short circuit



5. Double Fans  
Bring faster cooling and more efficiency sterilization





#### 6. Water Steam Separator

Patented design, prevent vacuum pump from oil damaging



#### 7. Brass Piping

Withstand high temperature and pressure safe and reliable



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)