



MEDICAL AUTOCLAVE BAMN-401

MEDICAL AUTOCLAVE BAMN-401

CLASS N STEAM STERILIZER

Designed for sterilization of solid instruments, medical class N autoclave is a non vacuum sterilizer which removes air from the chamber by thermodynamic air displacement. Safety features with a door locking system, these autoclaves are ideal for small spaces.



SPECIFICATIONS

Model	BAMN-401
Capacity	3 L
Sterilization Pressure	0.22 Mp
Sterilization Temperature	134 °C
Sterilization Time	4 min
Sterilization Cycle time	11 min
Required water	Distilled water or purified water
Product Size	300X190X350 mm
Outer Dimension	370x260x460 mm
Gross Weight	7.5 Kg
Power	0.5 KW
Voltage	AC220/110V+10%
Alt Name	Class N steam sterilizer

APPLICATIONS

Medical, Dental, Laboratory, Pharmaceutical, Research.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com