



MEDICAL AUTOCLAVE BAMN-403

MEDICAL AUTOCLAVE BAMN-403

CLASS N STEAM STERILIZER

Designed for sterilization of solid instruments, medical class N autoclave is a non vacuum sterilizer which removes air from the chamber by thermodynamic air displacement. Safety features with a door locking system, these autoclaves are ideal for small spaces.



SPECIFICATIONS

Model	BAMN-403
Capacity	18 L
Product Size	450X550X350 mm
Chamber Dimension	φ249x350 mm
Package Size	670x550x500 mm
Gross Weight	30 Kg
Net Weight	23 Kg
Power	800 W
Voltage	AC220/110V+10%, 50/60 Hz+ 1Hz
Alt Name	Class N Steam Sterilizer

APPLICATIONS

Medical, Dental, Laboratory, Pharmaceutical, Research.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com