



## MOISTURE ANALYZER BANA-504

# MOISTURE ANALYZER BANA-504

HALOGEN MOISTURE ANALYZER



5 Inch touch screen / CN and EN switch  
Highlight halogen lamp / HBM sensor  
Aluminum shell / Stainless steel chamber  
Real-time dynamic curve / Date and Time  
Temperature / Time set  
The percentage of moisture content / dry residual  
Store historical set  
Option: Aluminum box / Printer / Grinder / USB

## SPECIFICATIONS

Model	BANA-504
Capacity	110g
Readability	0.001g
Operating Temp	5°C-35°C
Display	5 Inch touch screen
Calibration	External calibration
Interface	RS232(Standard) USB(Optional)
Pan Size	Φ90mm
Packing Size	485*340*350mm
Gross Weight	6.8kg
Heating Unit	Halogen lamp-2
Temp.Sensor	PT100
Temp.Range	40°C-199°C
Temp.Step	1°C
Moisture Range	0.00%-100.00%
Moisture Readability	0.01%
Dry Range	100.00%-0.00%
Dry Readability	0.01%
Temp Setting	40°C-199°C By 1°C Step
Time Setting	1-99min By 10s
Number Of Storage	20
Number Of History	20
Heater power	220V ±15% 50Hz / 110V ±15% 60Hz / 400W
Balance power	Input:220V 50Hz / 110V 60Hz Output:9V
Alt Name	Halogen Moisture Analyzer

## ACCESSORIES FOR PURCHASE

No	Name	Size	Material	Mat
1	Aluminum sample tray 1	Outer circle Φ 90mm Inner circle Φ 83mm Plate height 4mm	Aluminum	50 pieces/box
2	Portable aluminum box	48*44*24cm		

3	Temperature kit		
---	-----------------	--	--



1



2



3

## FEATURES

- 5 Inch touch screen / CN and EN switch
- Highlight halogen lamp / HBM sensor
- Aluminum shell / Stainless steel chamber
- Real-time dynamic curve / Date and Time
- Temperature / Time set
- The percentage of moisture content / dry residual
- Store historical set
- Option: Aluminum box / Printer / Grinder / USB



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
 Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)