



## MOISTURE ANALYZER BANA-507

# MOISTURE ANALYZER BANA-507

INTELLIGENT MOISTURE ANALYZER



- 7 Inch touch screen
- Automatic door / Transparent window
- Highlight halogen lamp
- Polymer materials / Heat resistant molding
- Dual motor drive / Double fan design
- Temp. control / Stainless steel chamber
- Real-time dynamic curve / Date and Time
- Temperature / Time set
- Store historical set
- The percentage of moisture content
- Connect PC and printer
- Option: Printer / USB

## SPECIFICATIONS

Model	BANA-507
Capacity	210g
Readability	0.001g
Operating Temp	5°C-35°C
Pan Size	Φ100mm
Heating Source	Halogen lamp
Temp. Sensor	PT100
Temp. Range	40°C-230°C
Temp. Step	1°C
Moisture Range	0.00%-100.00%
Moisture Readability	0.01%
Dry Range	100.00%-0.00%
Dry Readability	0.01%
Temp Setting	40°C-230°C By 1°C Step
Time Setting	1-99min By 1s
Number Of Storage	20
Number Of History	20
Heating Unit	110V ±15% 60Hz/400W
Balance power	Input:220V 50Hz / 110V 60Hz Output:9V
Machine power	Input: 220V 50Hz / 110V 60Hz Output:24V Dual-drive
Host Dimension	415*215*195mm
Indicator Dimension	200*135*80mm
Operation Height	275mm
The Chamber Height	35mm
Display	7 Inch touch panel
Interface	RJ45/USB (optional)
Communication Mode	Cable
Rechargeable Battery	No
Alt Name	Intelligent Moisture Analyzer

## ACCESSORIES FOR PURCHASE

No	Name	Size	Material	Mat
1	Aluminum sample tray 2	Outer circle $\Phi$ 100mm Inner circle $\Phi$ 92mm Plate height 6mm	Aluminum	80 pieces/box
2	Temperature kit			



1



2

## FEATURES

- 7 Inch touch screen
- Automatic door / Transparent window
- Highlight halogen lamp
- Polymer materials / Heat resistant molding
- Dual motor drive / Double fan design
- Temp. control / Stainless steel chamber
- Real-time dynamic curve / Date and Time
- Temperature / Time set
- Store historical set
- The percentage of moisture content
- Connect PC and printer
- Option: Printer / USB



High brightness halogen lamp



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)