



# PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-104

# PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-104

BASIC LINE AUTOCLAVE



Durable & Reliable:

- a. Near 20 years of manufacture experience.
- b. Complete quality control system.


Effectiveness:

Stable control of temperature and pressure to guarantee the effectiveness during heavy and continuous sterilizations.

## SPECIFICATIONS

Model	BAPV-104
Capacity	23L
Input	220V±22V / 50HZ
Power	2300 VA
Chamber size	25.4*45 Cm
Product size	61*43.5*42.5 Cm(W*L*H)
Installation size	60Cm
N.W.	38Kgs
G.W.	40Kgs
Alt Name	Basic Line Autoclave

## FEATURES



Durable & Reliable:

- a. Near 20 years of manufacture experience.
- b. Complete quality control system.

Effectiveness:

Stable control of temperature and pressure to guarantee the effectiveness during heavy and continuous sterilizations.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)