



# PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-204

# PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-204

PRE AND POST VACUUM CLASS B AUTOCLAVE



- 1.Safety Valve  
Prevent chamber from over pressure
- 2.Bacteria Filter  
High quality,2µm precision
- 3.Steam Generator  
Patented design, high efficiency and reliable
- 4.Air Switch  
Shut down the machine in case of electricity leakage or short circuit
- 5.Double Fans  
Bring faster cooling and more efficiency sterilization
- 6.Water Steam Separator  
Patented design, prevent vacuum pump from oil damaging

## SPECIFICATIONS

Model	BAPV-204
Description / V/Hz	220/50
Voltage / W	2500
Capacity / L	22
Chamber size / cm	Φ 25x45
Product size / cm	66x50x43
Installation Size / cm	85
Weight / kg	55
Alt Name	Pre and Post Vacuum Class B Autoclave



## FEATURES



- 1.Safety Valve  
Prevent chamber from over pressure



- 2.Bacteria Filter  
High quality,2µm precision



3.Steam Generator  
Patented design, high efficiency and reliable



4.Air Switch  
Shut down the machine in case of electricity leakage or short circuit



5.Double Fans  
Bring faster cooling and more efficiency sterilization



6.Water Steam Separator  
Patented design, prevent vacuum pump from oil damaging



7.Brass Piping  
Withstand high temperature and pressure safe and reliable



8. German vacuum pump (optional)  
More powerful and less noise.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)