



# PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-501

# PRE AND POST VACUUM CLASS B AUTOCLAVE BAPV-501

## AUTOCLAVE (STEAM STERILIZER) CLASS B

Reliable and robust class B autoclaves intended to deliver the highest quality of sterilization. These autoclaves are easy to clean and its compact dimensions make it an ideal option for today's contemporary clinic.



Intelligent operating interface:

Intelligent, one-stop interface; and make the process of operation and sterilization clear at a glance.

Intelligent water cooling system:

Vacuum pump water cooling, making its long life strong and powerful.

Automatic door lock system:

It is easy to operate and completely eliminate the potential danger caused by improper operating of manual door closing.

## SPECIFICATIONS

Model	BAPV-501
Tank Capacity	45 L
Sterilization Temperature	121°C/110Kpa, 134°C/210Kpa
Vacuum degree	-096 Kpa
Residual humidity	<0.2%
121°C Solid	Yes
121°C Porous	Yes
121°C Hollow	Yes
134°C Solid	Yes
134°C Porous	Yes
134°C Hollow	Yes
B&D Test	Yes
Helix Test	Yes
Vacuum Test	Yes
Special sterilization	Inactivate AIDS (HV), hepatitis B (HBV) mad cow virus and spores
Drying program	Strong vacuum drying, the remaining humidity of the instrument <0.2%
Chamber Dimension	Ø319x620 mm
Overall Dimension (LxWxH)	845x605x540 mm
Packing Size	970x730x660 mm
Safety valve	Yes
Pneumatic lock	Yes
Fault self-check system	Yes
Sterilization record	Equipped with printer and USB interface to record the sterilization process
Display	LCD display, touch button
Door	Automatic locking mechanism
Weight	113 kg
Power	3500 W
Power Supply	AC220V/110V±10%
Printer	50/60Hz±1Hz
Alt Name	Autoclave (Steam Sterilizer) Class B



## FEATURES



Intelligent operating interface:

Intelligent, one-stop interface; and make the process of operation and sterilization clear at a glance.



Intelligent water cooling system:

Vacuum pump water cooling, making its long life strong and powerful.



Automatic door lock system:

It is easy to operate and completely eliminate the potential danger caused by improper operating of manual door closing.

## APPLICATIONS

Medical, Dental, Laboratory, Pharmaceutical, Research.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)